



# Cell SELEX

## Troubleshooting Guide

iGEM  
Stockholm  
2021





# INDEX

---

<b>BEFORE THE EXPERIMENTS</b>	<b>1</b>
DESIGNING AND ORDERING A SELEX LIBRARY	1
How to order a SELEX library on IDT - step by step	1
DESIGNING AND ORDERING PRIMERS	4
Forward primer	4
Reverse primer	4
How to order primers on IDT - step by step	4
OTHER NECESSARY PRODUCTS	8
<b>CELL SELEX PROTOCOL</b>	<b>10</b>
LIBRARY PREPARATION	10
NEGATIVE SELECTION	10
POSITIVE SELECTION	11
ELUTION OF BOUND APTAMERS	12
DNA PRECIPITATION	12
APTAMER LIBRARY AMPLIFICATION	13
SEPARATION OF SENSE AND ANTISENSE SSDNA	17
FLOW CYTOMETRY AND YIELD EVALUATION	18
<b>RELATED PROTOCOLS</b>	<b>19</b>
CFU DETERMINATION	19
2% AGAROSE GEL IN TBE	19
<b>BUFFERS AND SOLUTIONS</b>	<b>20</b>
BUFFERS	20
SOLUTIONS	21









## II. DESIGNING AND ORDERING PRIMERS

Primers are used to amplify by PCR the aptamers that bound to the target during the selection round.

Both primers need to include modifications in their 5' end in order to be able to differentiate between the sense and antisense DNA strand after the PCR and select the sense ssDNA strand.

### Forward primer

- Same sequence as the 5' overhang.
- 5' Fluorescent dye modification, i.e. 6-FAM (Fluorescein) → will label the sense strand.

I.e. our Forward primer sequence is: /56-FAM/ATC CAG AGT GAC GCA GCA

### Reverse primer

- Sequence is the reverse complementary sequence of the 3' overhang
- 5' Biotin modification → will label the antisense DNA strand.

I.e. our Reverse primer sequence is: /5BiosG/ACTAAGCCACCGTGTCCA

Taking advantage of the Streptavidin - Biotin affinity we can remove the antisense DNA strand from the library that will be used in the next SELEX round. I.e. by using Streptavidin beads.

### How to order primers on IDT - step by step

1. Enter on the IDT website (<https://eu.idtdna.com/pages>)
2. Click on Products and Services > Custom DNA/RNA oligos

The screenshot shows the IDT website header and a navigation menu. The IDT logo is on the left, followed by a search bar and utility links like 'GET HELP', 'EN', and 'STOCKHOLM IGEM 2021'. Below the search bar is a link to 'Order by stock part number'. The navigation menu includes 'PRODUCTS & SERVICES', 'APPLICATIONS & TECHNOLOGIES', 'SUPPORT & EDUCATION', 'TOOLS', and 'COMPANY'. Under 'PRODUCTS & SERVICES', there are four columns of links: 'COVID-19 SOLUTIONS', 'DNA & RNA' (circled in red), 'GENES & GENE FRAGMENTS', and 'NEXT GENERATION SEQUENCING'.

COVID-19 SOLUTIONS	DNA & RNA	GENES & GENE FRAGMENTS	NEXT GENERATION SEQUENCING
CDC assay	Custom DNA oligos	Double-stranded DNA fragments	Adapters
Charité/Berlin assay	Custom RNA oligos	Single-stranded DNA fragments	Library preparation
Luminex assay	Affinity Plus DNA & RNA oligos	Custom gene synthesis	Hybridization capture
Custom detection panels	DNA oligo pools	DNA Origami	Amplicon sequencing
NGS solutions	Large-scale synthesis		IDT Align Program
Cas13 guide RNAs	SameDay oligos		xGen Exome Research Panel v2
Affinity Plus ASOs	Inventoried oligos		
Genes & gene fragments	Oligo modifications		
Media information			



3. DNA oligos > Order now > Order in tubes

4. In Oligo Entry, increase the “# of items” to 2. Then, write the names and sequences of your primers and the scale (mol) of the product.

⚠ *Make sure that you are writing the Reverse primer sequence correctly: reverse complementary sequence of the 3' overhang, but also from 5' to 3' end.*

Purification > Standard Desalting

Select All    ACTIONS: ▾    # of Items: 2    GO    BULK INPUT 📄

# 1 Forward primer \* ⓘ

**Scale** ⓘ  
100 nmole DNA oligo ▾

**Sequence** \* (5' → 3')  
5' MOD ▾ INTERNAL ▾ 3' MOD ▾ BASES ▾  
ATC CAG AGT GAC GCA GCA

# Bases: 18 (Min:10 Max:90)    Min Yield: 35 nmoles  
GC: 55.6%    Tm: 56.5°C    ⚙ DeltaG: -34.45 kcal/mole

**Formulation**  
None ▾

**Purification**  
Standard Desalting ▾

**Services**  
🔍  
 Analytical RP-HPLC 420,00 kr  
 Analytical IE-HPLC pH 12.0 420,00 kr  
 Na+ Salt Exchange 570,00 kr

# 2 Reverse primer \* ⓘ

**Scale** ⓘ  
100 nmole DNA oligo ▾

**Sequence** \* (5' → 3')  
5' MOD ▾ INTERNAL ▾ 3' MOD ▾ BASES ▾  
|ACTAAGCCACCGTGCCA

# Bases: 18 (Min:10 Max:90)    Min Yield: 35 nmoles  
GC: 55.6%    Tm: 56.2°C    ⚙ DeltaG: -34.87 kcal/mole

**Formulation**  
None ▾

**Purification**  
Standard Desalting ▾

**Services**  
🔍  
 Analytical RP-HPLC 420,00 kr  
 Analytical IE-HPLC pH 12.0 420,00 kr  
 Na+ Salt Exchange 570,00 kr

5. Add the 5' modifications by clicking on "5' MOD"

- For the Forward primer, choose "/56-FAM/"
- For the Reverse primer, choose "/5BiosG/"

# 1 Forward primer

Scale 100 nmole DNA oligo

Sequence (5' → 3')

5' MOD INTERNAL 3' MOD BASES

/56-FAM/ATC CAG AGT GAC GCA GCA

# Bases: 18 (Min:10 Max:90) Min Yield: 35 nmoles

GC: 55.6% Tm: 56.5°C DeltaG: -34.45 kcal/mole

Formulation: None

Purification: Standard Desalting

Services

- Analytical RP-HPLC 420,00 kr
- Analytical IE-HPLC pH 12.0 420,00 kr
- Na+ Salt Exchange 570,00 kr

# 2 Reverse primer

Scale 100 nmole DNA oligo

Sequence (5' → 3')

5' MOD INTERNAL 3' MOD BASES

/5Biosg/ACTAAGCCACCGTGCCA

# Bases: 18 (Min:10 Max:90) Min Yield: 35 nmoles

GC: 55.6% Tm: 56.2°C DeltaG: -34.87 kcal/mole

Formulation: None

Purification: Standard Desalting

Services

- Analytical RP-HPLC 420,00 kr
- Analytical IE-HPLC pH 12.0 420,00 kr
- Na+ Salt Exchange 570,00 kr

6. Add to order (right side of the screen)

> Continue

Duplex | RxnReady | Plates

1 Items

ADD TO ORDER

ADD TO WISH LIST

Help

7. Double check:

- 1) Purification > Standard Desalting
- 2) Sequence: check if your sequence, modifications and sequence length is correct

# 1 Forward primer				ACTIONS ▾	🗑️	qty 1	GO	379,36 kr
<b>Product</b>	100 nmole DNA Oligo	<b>Expected Ship Date</b>	2021-10-06					
<b>Purification</b>	Standard Desalting	<b>Guaranteed Yield</b>	7.1 ODs = 35 nmol = 211.9 µgrams					
<b>Length</b>	18							
<b>Sequence</b>	/56-FAM/AT CCA GAG TGA CGC AGC A							

# 2 Reverse primer				ACTIONS ▾	🗑️	qty 1	GO	288,36 kr
<b>Product</b>	100 nmole DNA Oligo	<b>Expected Ship Date</b>	2021-10-06					
<b>Purification</b>	Standard Desalting	<b>Guaranteed Yield</b>	6 ODs = 35 nmol = 203.8 µgrams					
<b>Length</b>	18							
<b>Sequence</b>	/5Biosg/AC TAA GCC ACC GTG TCC A							

8. Check out.

Fill in your Shipping and Billing details and finalise the order.

CHECK OUT

ADD TO WISH LIST

E-MAIL CART / QUOTE

CONTINUE SHOPPING

### III. OTHER NECESSARY PRODUCTS

Besides the library and the primers, other products should also be order before starting with the SELEX experiments:

- Bacterial strains or Eukaryotic cell lines - Target of interest
  - ↳ for Positive selection
- Another bacterial strain from the same species / similar cell line
  - ↳ for Negative selection
- Nuclease-free water
  - ⚠ Nuclease-free water is used in multiple steps. In order to avoid contamination we recommend to prepare aliquots in 1.5 mL Eppendorf tubes and 15 mL Falcon tubes.
- Sodium acetate (NaOAc)
  - ↳ for DNA precipitation
- Ethanol 95%
  - ↳ for DNA precipitation
  - ⚠ Since it needs to be cold for the DNA precipitation, we recommend to store an aliquot in a 50mL Falcon at -20°C.
- Glycoblue
  - ↳ for pellet visualisation after DNA precipitation
- Taq polymerase
  - ⚠ In order to amplify and slightly modify the sequences in the library after the selection step in each SELEX round, a polymerase that lacks error correction (exonuclease 3'-5' activity) is needed.
- dNTPs
  - ↳ for PCR
- Agarose
  - ↳ for DNA electrophoresis
- DNA gel stain (I.e. SybrSafe)
  - ↳ for staining agarose gels in DNA electrophoresis
- Low range DNA ladder
  - ↳ for DNA electrophoresis

- Streptavidin beads (i.e. Dynabeads)
  - ↳ for separating the 5' 6-FAM sense strand from the 5' Biotin antisense strand after the PCR amplification
- Magnet rack
  - ↳ for separation of sense from antisense ssDNA bound to Dynabeads
- Sodium hydroxide (NaOH)
  - ↳ for eluting the sense ssDNA from the antisense ssDNA bound to the Dynabeads
- GE healthcare Cytiva illustra NAP column NAP-5
  - ↳ for desalting the library after eluting the sense ssDNA from the antisense ssDNA (bound to the Dynabeads) with NaOH
- MES free acid
  - ↳ for MES buffer
- Tris-base
  - ↳ for TBE buffer
- Boric acid
  - ↳ for TBE buffer
- EDTA
  - ↳ for TBE buffer
- PBS
  - ↳ for Buffer, solutions and washes
- BSA
  - ↳ for Wash Buffer
- MgCl<sub>2</sub>
  - ↳ for Wash Buffer
- Yeast tRNA
  - ↳ for Binding Buffer

**!** *Yeast tRNA can be purchased as a powder or liquid form. In case of acquiring the yeast tRNA powder, after resuspension, the solution should be frozen using liquid N<sub>2</sub> or dried ice and stored at -80°C.*


# CELL SELEX PROTOCOL


---

The SELEX protocol is very complex. Therefore, we went through every step and included a brief explanation of the rationale behind them as well as troubleshooting tips.


Although every type of cell (bacteria, yeast, eukaryotic cell) can be used to develop aptamers by cell SELEX, we will focus in bacterial Cell SELEX.

## LIBRARY PREPARATION

 Set the heater to 95°C in advance.


 Calculate how much Binding Buffer you need (~ 2.5 mL) and prepare 0.5 mL extra, just in case.

1. Let ssDNA pool denature by heating at 95°C for 5 min

 This step allows DNA denaturation and breaks down secondary structures.

2. Dissolve 20 pmol naive ssDNA library in 500 µl Binding buffer (BB)

▶ 20 µL of 1 µM library in 500 µL BB


 Make sure you don't need a big volume of you library. Otherwise it will dilute the BB.


3. Cool on ice for 10 min


 This step allows the formation of stable tertiary structure.


 Stop point

## NEGATIVE SELECTION → Don't do in the first SELEX round

 Prior to the negative selection, the cells needed for this step should be cultured. I.e. for bacteria, an O/N culture should be prepared the day before.


 CFU determination is required before starting the protocol. See Related protocols > CFU determination or look for OD600-CFU correlation in the literature (preferable).


 Very little volume of bacterial suspension is needed for selection steps. Prepare small cultures.


 *Don't expect to see proper pellets after the centrifugation in the washing steps. Remove the supernatant carefully and avoid touching the bottom and walls of the eppendorf.*


4. Wash  $\sim 10^7$  CFU "negative control" bacteria 2x 3 min 4000 rpm with 500  $\mu$ l Wash buffer (WB)
5. Wash bacteria 1x 3 min 4000 rpm with 500  $\mu$ l BB
6. Incubate ssDNA library with bacteria  $10^7$  CFU in 1mL BB for 30min on rotator at 4°C
7. Wash pellet 1 x 3 min 4000 rpm with 500 $\mu$ l WB
  - Supernatant contains unbound DNA sequences → New library **!**

## POSITIVE SELECTION


 *Prior to the positive selection, the cells needed for this step should be cultured. I.e. for bacteria, an O/N culture should be prepared the day before.*

 *CFU determination is required before starting the protocol. See Related protocols > CFU determination or look for OD600-CFU correlation in the literature (preferable).*


 *Very little volume of bacterial suspension is needed for selection steps. Prepare small cultures.*

 *Don't expect to see proper pellets after the centrifugation in the washing steps. Remove the supernatant carefully and avoid touching the bottom and walls of the eppendorf.*

8. Wash  $\sim 10^7$  CFU "positive" bacteria 2x 3 min 4000 rpm with 500  $\mu$ l Wash buffer (WB)
9. Wash "positive" bacteria 1x 1min 4000 rpm in 500  $\mu$ l BB
10. Incubate ssDNA library (supernatant with unbound sequences) with  $10^7$  CFU bacteria in 1mL BB for 1h on rotator at 4°C
11. Wash pellet 1 x 3 min 4000 rpm with 500 $\mu$ l WB

 *In later rounds: increase number of washing steps by one more wash and one more minute, until you reached 5 minutes of washing time.*


## ELUTION OF BOUND APTAMERS

 Set the heater to 95°C in advance.

12. Resuspend in 500µL BB
13. Heat at 95°C for 15min
14. Cool on ice for 5 min
15. Centrifuge at 14,000 rpm for 2 min → Supernatant contains aptamers

 Freeze point

## DNA PRECIPITATION

 Remember to have cold 95% EtOH → -20°C

16. Add 0.1V x 3M Sodium acetate to the supernatant from step 19
  - ▶ 50 µL 3M Sodium acetate for 500 µL
17. Add 100µg/mL Glycoblue to solution
  - ▶ 5 µL Glycoblue for 500 µL
18. Add 2.5V 95 – 99% EtOH, cold (- 20 °C)
  - ▶ 1.375 mL EtOH for 500 µL
  - a) Incubate 1h at -20°C
  - b) Incubate O/N at -80°C
19. Centrifuge at 10.000g for 30 min at 4°C
20. Remove supernatant
21. Wash in 50 µL 70% EtOH
  - ⚠ The EtOH volume should be the just enough to cover the pellet.
22. Take as much supernatant as possible and let dry
  - ⚠ You shouldn't be able to smell EtOH anymore. Shouldn't look glossy. Don't let it dry too much either otherwise pellet resuspension becomes harder.
23. Resuspend in 30 µL of nuclease free water
24. Measure DNA concentration with Nanodrop

 Freeze point



## APTAMER LIBRARY AMPLIFICATION

### I. Library pre-amplification

29. Prepare 1 × 50  $\mu$ L PCR tube:

Component	50 $\mu$ L reaction	Work concentration
Forward primer (10 $\mu$ M)	5 $\mu$ L	1 $\mu$ M / 50 pmol
Reverse primer (10 $\mu$ M)	5 $\mu$ L	1 $\mu$ M / 50 pmol
dNTP mix (2.5 mM each)	5 $\mu$ L	0.25 mM each
DNA template (~100 ng/ $\mu$ L)	1 $\mu$ L	100 ng
Taq DNA Polymerase	-	1.25 U
Polymerase Buffer	-	1X
H2O (nuclease-free)	Up to 50 $\mu$ L	-

⚠ *Always keep tubes on ice*

⚠ *Add polymerase last*

⚠ *Check your polymerase product sheet to calculate how much DNA pol and DNA pol buffer are needed.*

⚠ *Use filtered tips and be careful during the pipetting to avoid contamination.*

30. PCR protocol for Library pre-amplification:

Initial denaturation → 95°C 3'

-----

Denaturation → 95°C 30''

Annealing → -°C 30'' × 6 Cycles

Elongation → 72°C 30''

-----

Final elongation → 72°C 3'

⚠ Annealing temperature → General rule: the annealing temperature should be 3–5°C lower than the lowest T<sub>m</sub> of the primers.

⚠ Elongation time → For Taq polymerase, elongation time is 1 min per 1 kb DNA, minimum elongation time being 30''.

⚠ Temperatures and times in the PCR protocol were selected based on our Taq polymerase. Always check the recommended protocol in the User manual from your manufacturer and adapt them accordingly.

! Now you have a new library!!!

❄ Freeze point

## II. Library amplification

31. Prepare 4 tubes × 50 µL PCR tube:

- ▶ 3 tubes for the cycle testing → optimising the amount of cycles for your library is needed. Each tube will be run for different amount of cycles and the resulting amplifications will be compared.
- ▶ 1 tube for NTC (Non Template Control) → 0 µL DNA template

Component	50µL reaction	Work concentration
Forward primer (10µM)	5 µL	1 µM / 50 pmol
Reverse primer (10µM)	5 µL	1 µM / 50 pmol
dNTP mix (2.5 mM each)	5 µL	0.25 mM each
DNA template (~100 ng/µL)	1 µL	100 ng
Taq DNA Polymerase	-	1.25 U
Polymerase Buffer	-	1X
H2O (nuclease-free)	Up to 50 µL	-

⚠ *Always keep tubes on ice.*

⚠ *Add polymerase last.*

⚠ *Use filtered tips and be careful during the pipetting to avoid contamination.*

⚠ *Preparing a Master mix also prevents contamination and is less time consuming.*

32. PCR for Library amplification:

Initial denaturation → 95°C 3'

-----

Denaturation → 95°C 30''

Annealing → -°C 30'' × 35 Cycles

Elongation → 72°C 30''

-----

Final elongation → 72°C 3'

- ▶ Take out one tube after: 20, 25, 30 and 35 cycles

! **Annealing temperature** → General rule: 3–5°C lower than the lowest T<sub>m</sub> of the primers.

! **Elongation time** → For Taq polymerase, elongation time is 1 min per 1 kb DNA, minimum elongation time being 30''.

! Temperatures and times in the PCR protocol were selected based on our Taq polymerase. Always check the recommended protocol in the User manual from your manufacturer and adapt them accordingly.

❄ *Freeze point*

33. Run 2 μL of each PCR product on a 2% agarose, 120V for 15 min

! See 2% agarose gel in TBE

- ▶ 2% agarose gel: 1 g agarose, 50 mL 1% TBE, 5 μL SybrSafe
- ▶ Low range DNA ladder: 2 μL Ladder (ready to use, dye not needed), 4 μL ddH<sub>2</sub>O
- ▶ Samples: 2 μL sample, - μL DNA dye (working concentration 1X), - μL ddH<sub>2</sub>O (up to 6 μL)

↪ Decide the **optimal PCR cycle** number based on the smearing in the gel. Choose number of cycles from the sample with less smeared and more clear band.

34. Repeat PCR library amplification (Protocol step 31, 32) with the number of cycles determined in step 33:

6 cycles (from pre-amplification) + ??? cycles (from amplification)

⚠ Do not prepare only one PCR tubes. Alternatively, prepare a Master Mix for several reactions in order to have a big amount of library, which is needed for the next step. It will also be helpful in case troubleshooting is needed and because.


### ❄ Freeze point

📌 ! Based on our experience, optimising the PCRs in this step could take some time. Here we gathered some **extra tips** that could save you time along this process:

- We recommend spinning down the tubes before placing them in the thermocycler, to make sure that every drop is in the mix.
- Be extremely careful in order to avoid contamination of any of the PCR reagents.

**If contamination appears**, prepare a new PCR tubes were you change one component at a time: i.e. tube 1 - new dNTPs, tube 2 - new Fwd primer, etc. Then run a gel and interpret the bands: the undesired band(s) will disappear in the sample where you replaced the contaminated component.

- **If there is no (or little) product in the gel and you can only see primer dimers:**





On the left image, all the bands correspond to primer dimers. On the right image, clear bands that correspond to the right size of our aptamer (~ 70 bp)

- 1) Check the final concentrations of every component. Your stocks might differ from ours, thus you cannot use the volumes we suggest in the tables. Look at the “Work concentration” column and adjust your volumes.
- 2) Check the annealing temperature. When primers are not annealing properly to the template DNA the PCR product will consist of primer dimers and little to no desired product. Check the annealing temperature of your primers in their respective IDT sheets. The annealing temperature should be 3–5°C lower than the lowest T<sub>m</sub> of the primers. Using T<sub>m</sub> calculator tools might be helpful, i.e. [www.thermofisher.com/tmcalculator](http://www.thermofisher.com/tmcalculator)

- We recommend **optimising the PCR in pre-amplification step**. There is no need to waste reagents and time with following PCRs, you can continue once the PCR is working properly.
- If the optimisation is not being straightforward, **don't waste your eluted library** (previously incubated with the targets). Use the raw library instead and once the PCR is working move to the eluted library to continue with the SELEX cycle.
- **If you don't see clear bands** (there's still smear) when running "20, 25, 30 and 35 cycles" samples in the gel, check lower and higher number cycles. The optimal one might be before 20 or after 35!

## SEPARATION OF SENSE AND ANTISENSE SSDNA

 *In this step using the right proportion beads - DNA is key.*

 *Measuring the DNA concentration with Nanodrop of a PCR product is not reliable. The sample contains PCR byproducts, primers, dNTPs, polymerase etc. In order to estimate the real DNA concentration in your sample you should compare the bands from the ladder and the PCR product in a gel image. Since the DNA concentration and loaded volume in the ladder are known, it is possible to make a rough DNA concentration estimation.*

3. Add the Streptavidin beads to the DNA


 *Prepare, wash and incubate the DNA with the beads according to the manufacturer.*

4. Add the magnets and discard the supernatant

5. Wash the beads + bound sequences with 500  $\mu$ L WB

6. Elute the ssDNA from the beads by melting in a 0.1 M NaOH solution during 5'

7. Add the magnets and recover supernatant

 *Try to keep the volume under 600  $\mu$ l. Keep in mind this will be your library for the next round and you don't want the aptamers to be very diluted.*

8. Desalt by using GE healthcare Cytiva illustra NAP column NAP-5

- Run the protocol as given in the kit to filter out the NaOH
- Recover the flow through

## FLOW CYTOMETRY AND YIELD EVALUATION

9. Repeat binding step (see *Positive selection*)
10. Flow cytometry to assess library enrichment after the SELEX round
  - ⚠ *There is no need to perform flow cytometry in every SELEX round, especially at the beginning, since the increase in affinity will not be that big. It is recommended to do it every 2-3 rounds.*

📌 SELEX cycle 1 finished!


Repeat until the flow cytometry results show a plateau in the fluorescent intensity. Then, send the “plateau library” for sequencing to know your aptamers sequence. Probably, there will be different DNA species. Next step would be to perform a binding assay with each one of them separately to determine their  $K_d$  and choose the one with higher affinity.

# RELATED PROTOCOLS

---

In this section protocols related to the SELEX cycle are included.

## CFU DETERMINATION

1. Take a 1 mL sample of an O/N bacterial culture
2. Create a dilution series: at least 10 series of 1:10 dilution
3. Plate 100  $\mu$ L of each series onto LB Agar plates (+ Antibiotic)  
 *Media and antibiotic need to be chosen according to your bacteria*
4. Incubate the plates overnight in the culture conditions indicated for your strain
5. Select the plate with a countable amount of colonies
6. Count the colonies and calculate the CFU/mL:
  - ▶ Each colony corresponds to 1 CFU
  - ▶ Calculate the CFU/mL taking into account the inoculation volume (100  $\mu$ L) and the specific dilution factor of your plate (1:10, 1:100, 1:1000, ...)

## 2% AGAROSE GEL IN TBE

1. Prepare a 2% agarose solution in TBE in a beaker (or similar glass recipient).
2. Heat the beaker in the microwave until the solution boils.
3. Stir and heat again until reaching the boiling point again.
4. Stir and make sure all the agarose is dissolved and the solution is clear.
5. Add the DNA gel stain.
6. Stir until it is dissolved homogeneously.
7. Pour the solution into the cast (small 50 mL cast) and put the comb.
8. Remove bubbles with a tip and wait ~30 min for the gel to polymerise.

# BUFFERS AND SOLUTIONS

---

We recommend to filter all buffers and solutions using a 0.2  $\mu\text{m}$  filter to prevent contamination and maintain the samples' purity.

## BUFFERS

### Wash Buffer (WB)

- PBS (1X)
- BSA (1 mg/mL)
- 5mM  $\text{MgCl}_2$

### Binding Buffer (BB)

- yeast tRNA (0.1mg/mL)  
in WB

### MES buffer (0.5M pH=6)

- MES free acid
- $\text{dH}_2\text{O}$

### TBE buffer (pH=8.5)

- Tris-base 0.89 M
- Boric acid 0.89 M
- EDTA 25 mM
- $\text{ddH}_2\text{O}$



## SOLUTIONS

### Sodium acetate (3 M)

- NaOAc
- ddH<sub>2</sub>O

### Sodium hydroxide (0.1 M)

- NaOH
- ddH<sub>2</sub>O